

1/3

Figure 1





Figure 2

## Detection of A Fluorescent Signal On A Lateral Flow Test Strip



**Figure 3 — Table of Exemplary Fluorescence Excitation Sources**

Source	Principal Lines (nm)
Mercury-arc lamp	366, 405, 436, 546, 578 *
Xenon-arc lamp	250–1000 *

Tungsten-halogen lamp	350–1000 *
Cr:LiSAF laser	430
Helium-cadmium laser	325, 442

Argon-ion laser	488, 514
Nd:YAG laser	532
Helium-neon laser	543, 594, 633

Krypton-ion laser	568, 647
* Continuous white-light source. Frequency-doubled principal line output.	

---